



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ VARIMED Treatment table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2018,)
- Upholstery with medical approval in orange color (upholstery can be replaced with a new one from our color chart for an additional fee),
- SCHMITZ VARIMED treatment table has a wide range of height adjustment for ergonomic work at the right height. The specially shaped lower part provides sufficient legroom,
- Technical specifications:
 - Manual adjustment:
 - Height range: 49 - 91 cm,
 - Headrest inclination: -33° / 0° / +33°,
 - Driving wheel guidance system,
- Mobile treatment table on 4 braking casters (built-in central brake),
- Includes:
 - Removable paper rack,
 - Upholstery has normal signs of use,
- Dimensions:
 - Lying surface: 66 x 144 x 8 cm,
 - Headrest: 66 x 51 x 8 cm,
 - Overall: 66 x 197 x 49/91 cm,
- Maximum load capacity: 225 kg,
- Weight: 98 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

