



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS INTELLIVUE MP30 Patient monitor / cardiac monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 1 built-in 11.1 V / 7.8 Ah battery (can be retrofitted with a second one at additional cost),
- Frequency: 50/60 Hz,
- IntelliVue MP30 is a cardiac monitor that monitors multiple vital functions,
- It has:
 - Panel screen, made with TFT LCD technology with a diagonal of 10" and working in SVGA resolution,
 - Touch screen operation, and navigation joystick, this allows direct access to many functions with simple commands,
 - Convenient handle and sturdy case make it easy to carry,
 - Trends of measurements can be presented graphically, combining up to three curves on a single graph, making it easy to identify changes in the patient's physiological state,
 - Trends of fast-variable parameters can be updated every single beat and up to four high-resolution trend segments can be displayed,
 - IntelliVue technology - access information from various systems in the hospital network - pharmacy, laboratory, ECG management, radiology department PACS, server-based hospital systems and others - and deliver it to the patient monitor,
 - Lacking a separate hard drive and fan, IntelliVue monitors are small, quiet, stable and quick to boot up,
 - Navigation dial for easy entry A navigation dial makes it easy to enter information and navigate the screen,
 - Audible and visual alarms,
 - LAN,
- Equipped with PHILIP Model M3001A module with wiring:
 - ECG /RESP - 5 - lead cable with extension,
 - SPO2 - Brand new sensor,
 - NBP - New cuff with air tubing,
 - TEMP - Temperature connector socket, (no sensor included, can be retrofitted at additional cost),
- Menu in language: Polish,
- Included:
 - English user manual (PDF),
 - 5-lead ECG cable,
 - Brand new NIBP cuff,
 - Brand new SPO2 sensor,
 - Power cord,
- Dimensions: 21 x 34 x 28 cm,
- Weight: 6.5 KG,
- Cat 02,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

