



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SILFRADENT APAG Surgical syringe heater (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Italy, (2017),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 100 VA,
- APAG is a medical device for heating syringes at a controlled temperature depending on the program,
- The goal is to obtain a biologically compatible polymeric material used as a filler material and to stimulate connective tissue to produce collagen used in aesthetic medicine,
- The secondary and tertiary structures are denatured, while the peptide bonds of the primary structure between amino acids remain intact,
- The density of the product depends on the heating time and is adjusted to the site of administration. The specially designed medical device is anti-static/antimagnetic. It is safe in the process of heating the product,
- Technical specifications:
 - Temperature control: 35°C to 80°C,
 - Time adjustment: 1 min. - 60 min,
 - 2 preset programs:
 - Albumin denaturation (time adjustable),
 - Collagen activation (time adjustable),
 - Built-in 3 programs with freely adjustable temperature and time,
 - Built-in heating cups:
 - 4 pcs. - ø11 mm. x 110 mm,
 - 4 pcs. - ø8.4 mm. x 110 mm,
- Selectable menu in language: ENG, FR, GER, IT,
- Keypad on control panel in language: ENG,
- Included:
 - English user manual (PDF),
 - Power cord,
- Dimensions: 14 x 25 x 17 cm,
- Weight: 2.6 KG,
- Has an up-to-date inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

