



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STORZ 250 TWIN 20113320 Light source with endoscope set (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 330 VA,
- Storz 250 TWIN is a halogen light source suitable for flexible endoscopes,
- Technical specifications:
 - Halogen light,
 - 2 bulbs 24 V 250 W,
 - Color temperature: 3000 K,
 - Bulb life: 50 h,
 - 3-stage adjustment of light intensity with a knob,
 - Bulbs protected by additional cover,
 - Manual selection of the bulb to work with the included allen key hidden under the cover of the device,
 - Input for the other end of the fiber optic cable protects it from damage during a break in operation,
 - Screw-out fiber optic input allows its installation of an adapter corresponding to another manufacturer's fiber optics (adapter not included),
 - Foiled control panel,
- Cysto - rigid urethroscope KARL STORZ 27026CB:
 - Color code: red,
- Cysto - rigid urethroscope KARL STORZ 27026CB;
 - Color code: white,
- KARL STORZ 27072F Tissue grasping forceps,
- KARL STORZ 27072A tissue grasping forceps,
- 2 pcs. - KARL STORZ 27025F telescoping bridge,
- KARL STORZ 27182A flexible injection cannula,
 - Length: 50 cm,
 - 6 FR. - (2 mm.),
- KARL STORZ 27005 BA HOPKINS II rigid endoscope:
 - Working length: 30 cm,
 - Viewing direction: 30° anterior oblique telescope,
 - Color code: red,
 - Autoclavable,
- KARL STORZ cystoscope cover,
- Included:
 - Power cord,
 - 2 pcs. - STORZ 495 NA fiber optic cable,
 - Flexible monopolar coagulation electrode 7Fr. / 53 cm,
 - KARL STORZ 27025G telescopic bridge,
 - 2 pcs. - plugs,
 - 2 pcs. - Steel faucet with Luer Lock connection KARL STORZ 27502,
 - 2 pcs. - Luer Lock connector 27500,
 - 2 pcs. - Luer Lock connector 600007,
 - 2 pcs. - spare 15 V / 150 W bulbs,
 - KARL STORZ valve grease,
- Dimensions: 30.5 x 11 x 24.5 cm,
- Weight: 7.1 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

