



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OEV 262-H Medical monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan, (2017),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 112 VA,
- The OEV-262H high-resolution MONITOR is a 26-inch LCD monitor with LED backlight for optimal use of the Olympus endoscope system,
- Features:
 - Full HD 26-inch LCD panel with LED backlight, 16:9 aspect ratio, high brightness, high contrast and high gradation images,
 - FLIP functions (such as mirroring and 180° rotation) that provide useful display patterns and monitor layouts for procedures,
 - Multi-modal display capabilities, including Picture in Picture (PIP), Picture out Picture (POP) and Clone Out1 with various image size combinations,
 - AIME (Advanced Image Multiple Enhancer) improves image quality by increasing resolution and reducing noise, enhancing the viewing experience of the Olympus endoscopy system,
- Technical specifications:
 - 26" (inch) TFT-Si LCD matrix,
 - Resolution: 1,920 × 1080 pixels (full HD),
 - Contrast ratio: 1,400: 1,
 - Aspect ratio: 16: 9,
 - Viewing angle: 178 ° (horizontal and vertical),
 - Number of colors: 1.07 billion,
 - Input - output connectors:
 - VIDEO BNC: IN/OUT,
 - S-VIDEO Y/C: IN/OUT,
 - HD15/VGA D-sub 15-pin: IN/OUT,
 - DVI connector (1): IN/OUT,
 - DVI connector (2): IN/OUT,
 - VGA 9-pin connector: RS 232C,
 - BNC/SDI-1 connector: IN/OUT,
 - BNC/SDI-2 connector: IN/OUT,
 - BNC Aux in: (HD / SD-SDI) connector,
- Selectable language menu: ENG, CHIN, IT, ESP, GER, FR,
- Included:
 - Power cord,
- Dimensions: 62.8 x 79.6 x 39.5 cm,
- Weight: 7.8 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

