



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **RECK MOTOMED LETTO 2 Physical therapy device (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2014),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- MOTOMed letto 2 is equipped with large, easy to feel buttons so that the patient can reduce or increase the resistance when pedaling, adjust the speed of pedal movement,
- Application:

◦ Passive training using motor power - this type of exercise is recommended for spastic paralysis, as a prelude to further physiotherapy.

With passive exercises, the muscles of the immobilized patient relax and muscle tension is gradually reduced,

◦ Active training using motor and pedal assist - recommended for patients who can perform moderate exercises on their own. The motor and pedal-assist software make it easier for the user to perform the exercises correctly,

◦ Active training - based solely on the patient's muscle strength, with the ability to adjust different levels of difficulty (deceleration force).

Before and after active training, it is recommended to perform passive training as a warm-up and muscle relaxation,

- Features:
  - Push-button operated device,
  - Even pedaling thanks to the relaxation drive,
  - Built-in program for relaxing spastic conditions,
  - Automatic spasticity recognition function,
  - Symmetrical training function,
  - Measurements of muscle tension at different stages of training,
  - Training analysis:
    - Speed (rpm),
    - Degree of difficulty (muscle tension),
    - Distance covered (km),
    - Energy consumed (Joule/kcal),
    - Brake resistance (gear/grade of load),
    - Achievement data (Watt),
- Technical specifications:
  - Sturdy metal construction with expandable base,
  - Number of revolutions per minute adjustable from 0 to 60,
  - Control with illuminated remote control,
  - Transport rollers for easy movement,
  - Two-stage change of pedal sweep size,
  - Spacing between feet 16 cm,
- Selectable menu in ENG, FR, ESP, IT, PT, NED, DAN, NOR, PL, RUS, TURK, GER, ICE,
- Included:
  - User manual in Polish (PDF),
  - Power cord,
- Dimensions: 50 x 118 x 124 cm,
- Weight: 65 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

