



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GE MAC 1600 ECG unit (Ex-Display)

- New (ex-display),
- Real photos of the product,
- Made in Germany, (2014),
- Power supply: 230 V,
- Emergency power supply: via 2 built-in Li-Ion 14.8 V / 2.3 Ah batteries,
- Frequency: 50/60 Hz,
- Technical specifications:
 - Built-in clear "7" inch color display,
 - Compact housing,
 - Device can be configured to read resting ECG, rhythm, arrhythmia and stress tests,
 - Detects improper patient connection,
 - Built-in printer,
 - Recording speed: 5 / 12.5 / 25 or 50 mm/sec,
 - Recording sensitivity: 2.5 / 5 / 10 / 20 / 40 mm/mV,
- Selectable menu in EN, DE, FR, IT, ES, SE, PT, DU, FI, HU, PL, CS, NR, CN, JP, RU, SK, KR,
- Apparatus placed on a mobile cart from GE Healthcare:
 - 4 stable casters including 2 braking casters,
- Included:
 - Power cord,
 - New 10-lead ECG cable,
 - 2 x cable clips,
 - Mactrode 10 wide clip set,
 - Attachable ECG lead holder,
 - Arm with holder for attaching ECG leads,
- Total dimensions: 46 x 60 x 120 cm,
- Total weight: 29 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

