



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**ELLMAN PELLEVE S5 Radiofrequency surgical diathermy (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 120 W,
- Pelleve S5 is a device using two types of radio waves for precision surgery and skin lifting without recovery period,
- Features:
  - Pelleve precisely delivers energy to skin tissue using advanced radiofrequency technology to induce collagen contraction without damaging the epidermis,
  - Compact ergonomic design allows convenient maneuvering around facial contours, including hard-to-reach areas such as the periorbital area,
  - With precise cutting, wounds heal much faster and are hardly visible,
  - Can be used for various types of procedures,
  - Low risk of swelling after procedures,
  - Precise cutting in various tissues,
  - Digital displays of the set power,
  - Foiled, easy-to-use control panel,
- Uses two types of radio waves for precision surgery and skin lifting without recovery period:
  - 4.0 Mhz monopolar radio wave,
  - Bipolar radio wave with a frequency of 1.7 Mhz,
- It has the following functions:
  - Monopolar:
    - CUT - 120 W,
    - BLEND - 90 W,
    - COAG - 60 W,
    - FULGURATE - 45 WATT,
  - Bipolar:
    - CUT - 120 W,
    - BLEND - 120 W,
    - COAG - 120 WATT,
    - FULGURATE - 120 W,
- Diathermy placed on a mobile ELLMAN cart equipped with:
  - Convenient carrying handle,
  - Movable arm with clip,
  - Three shelves for accessories,
  - Four double casters, including two front casters with brakes,
- Includes:
  - Handles:
    - Pelleve GlideSafe 15mm. - 2 pcs,
    - Pelleve GlideSafe 10mm. - 3 pcs,
    - Ellman with switch - 3 buttons,
    - Ellman handle,
  - BALLET needle electrode,
  - Ellman 7.25" x 1.0 mm tweezers,
  - Ellman tweezers 3.5" x 0.5 mm.,,
  - ELLMAN flat electrode,
  - Foot switch,
  - Electrode set - 4 pcs,
  - Connection wires 3 - pin - 4 pcs,
  - Connecting cables 2 - pin - 2 pcs,
  - Power cord,
- Overall dimensions: 50 x 35 x 130 cm,
- Weight: 39.5 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

