



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDIKA Carboxytherapy machine (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland, (2016),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- CARBOXYTHERAPY MEDIKA device is designed for carbon dioxide therapy, which is a major breakthrough in aesthetic medicine due to its wide range of applications and high effectiveness. The treatment involves intradermal or subcutaneous injection of appropriate doses of carbon dioxide,
- Application:
  - Dark circles under the eyes,
  - Biorevitalization of the face,
  - Biorevitalization of the décolletage and neck,
  - Loose skin (arms, abdomen, thighs, buttocks),
  - Cellulite (arms, abdomen, thighs, buttocks),
  - Excessive fatty tissue (arms, abdomen, thighs, buttocks),
  - Stretch marks,
  - Psoriasis,
- Technical specifications:
  - Operating modes:
    - 23 ready-made treatment programs,
    - Manual mode - customizable parameters,
  - Touch screen: 8" inches,
  - Gas temperature control: 20°C - 50°C,
  - Gas flow control: 40cc/m To 200cc/m,
  - Gas dosage adjustment: From 5cc To 5000cc
- Systems:
  - GPC system - Gas pressure control - electronic gas pressure control system and use of medical CO2 regulator (maintenance-free regulator,
  - System PGH - Precise gas heating - gas heating system adjustable from 20°C to 50°C, significantly reduces pain sensation during the performed procedure,
  - AFC - Automatic flow control system - the automatic gas flow control system allows precise regulation of carbon dioxide dosage in the range of 40 cc to 300 cc,
  - AAV system - Automatic air venting - automatic air venting system of the device - complete elimination of air from the device system minimizes the risk of postoperative complications,
- Included:
  - Power cord,
  - Reducer,
  - CO2 cylinder,
- Dimensions: 40 x 40 x 97 cm,
- Weight of the device: 20 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

