



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



K2-230C K2 TABLES Ophthalmic table (Ex-Display)

- New (ex-display),
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2450VA,
- The K2 ophthalmoscope table provides a fixed and sturdy design that gives confidence and a sense of comfort. It is height adjustable - ACP technology ensures the highest standards of performance and reliability,
- Built-in 4 power cord outlets for free connection of devices on the table,
- Built-in plastic shelf for ophthalmic accessories with two pull-out drawers,
- Free electric height adjustment in the range: 79 - 100 cm,
- Mobile table on 4 casters (2 braking),
- Includes:
 - User manual in English,
 - 1 x new 110 V power cord,
 - 1 x new 230 V power cord,
 - 2 x new power cord instrumentation,
 - 4 x mounted feet for table base,
 - Set of plastic fasteners for wiring,
- Dimensions:
 - Tabletop: 91 x 61 cm,
 - Overall 91 x 61 x 94/115 cm,
- Maximum lifting capacity: 125 KG,
- Weight: 43 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

