



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**EMED SPECTRUM ARGO PLASMA Surgical diathermy (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Poland, (2020),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1350 VA,
- The EMED ARGO PLASMA is an electrosurgical coagulator that uses high-frequency current energy to produce a low-temperature argon plasma affecting tissue,
- One of the key features of argon plasma technology is the minimal risk of carbonisation during soft tissue coagulation. This makes procedures not only safe, but also easier to control. In this way, specialists can achieve a high degree of precision and efficiency in their procedures,
- Advantages:
  - The energy allows the coagulator to selectively act on the least resistant tissues - other structures remain undamaged,
  - The use of the coagulator ensures the safety and efficiency of the surgical process,
  - High level of precision and thermal control,
  - Uninterrupted operation even with prolonged and frequent switching,
  - Argon properties ensure effective cooling of the operating field,
  - Minimises the risk of fabric carbonisation and smoke formation,
  - Consistent minimum depth of thermal exposure - limited depth of tissue damage,
  - Accelerated regeneration of damaged tissue and shorter post-operative recovery period,
- Features:
  - Universal use of the coagulator - the device has all the functions of an electrosurgical device,
  - Possibility to use the coagulator for other surgical procedures,
  - Functional reusable control handle with buttons to activate RF current,
  - Ceramic electrode tip minimises heating of the instrument and facilitates access to tissue layers,
  - Argon plasma coagulation, cutting, monopolar or bipolar coagulation, ability to change and create modes/programs,
- Application;
  - Plastic surgery:
    - Liposuction procedures,
    - Bariatric surgery,
    - Skin tightening with rapid weight loss,
    - Acceptable during pregnancy,
  - Mammoplasty:
    - Reconstructive surgery in mammoplasty,
    - General surgery:
    - Widely used in laparoscopy and endoscopy,
- Coagulator modes:
  - Argon plasma coagulation,
  - Cutting,
  - Monopolar or bipolar coagulation,
- Technical specifications:
  - Monopolar CUT: 400 W,
  - Monopolar COAG: 200 W,
  - Bipolar CUT: 400 W,
  - Bipolar: COAG 350 W,
- Equipment:
  - EMED 2-bolt 29 mm bipolar cable,
  - Argon electrode holder, wide, 2 buttons, 3.5m cable, SDS plug EMED,
  - 4mm electrode holder, wide, 2 buttons, 4m cable, 3-pin EMED plug,
  - Bipolar cable for ThermoStapler instruments 3 m long with BISSINGER forceps,
  - Neutral silicone electrode with application tape, reusable, for adults 4m EMED,
  - Bipolar tweezers straight 30 cm x 2 mm,
  - Bipolar tweezers straight 23 cm x 2 mm,
  - 2 pcs. - Bipolar tweezers 22 cm x 2 mm,
- Included:
  - User manual in PL, EN, IT, PT, ES (PDF),
  - Power cord,
  - 7 x reusable rigid argon electrode for k
  - 7 x reusable rigid argon electrode for coagulation, working length 100 mm., 5 mm. diameter,
  - 29 x multi-use rigid argon electrode for coagulation, working length 320 mm., diameter 5 mm.,
  - The last photo shows the accessories together as a set,
- Dimensions: 58 x 56 x 130 cm,
- Weight: 44 KG,
- Has an up-to-date inspection and is ready for use,
- Technical Passport (Service Report) issued, valid for a period of 12 months,
- Warranty:
  - 6-months for the domestic market (Poland),
  - 3 months for the international market,
  - Possibility to extend the warranty to 12 months for the domestic and international market for a surcharge,
- Financing options (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

