



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET RADIUS 1557.02FJ Gynecological chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 290 VA,
- Maquet Radius is used for gynecological and urological examinations and procedures,
- Upholstery with medical approval in blue color (it is possible to replace the upholstery with a new one from our color pattern at an additional cost),
- Possibility of electromechanical adjustment by wired hand and foot remote control:
 - Bed height in the range of 61 - 91 cm,
 - Tilt angle:
 - Back segment in the range (- 9°) / 0° / 57°,
 - Seat in the range of 18°,
 - When adjusting the angle of the back segment, the seat angle is automatically adjusted,
 - Storage compartment for a 38 cm wide roll of backing paper,
 - Cutout in the upholstery of the back segment for the backing paper,
 - 3 pre-programmed buttons on the handheld remote control,
 - S button to stop the programmed button mode,
 - Pull-out waste bowl with Ø 21 cm,
 - Ability to electrically unfold the seat into a reclining position,
 - Ability to freely adjust the foot rest vertically and 45° horizontally,
- Includes:
 - User manual in Polish (PDF),
 - Power cord,
 - Wired handheld remote control,
 - Wired foot remote control,
 - Waste bowl,
- Dimensions:
 - Back segment: 91 x 61 cm,
 - Seat: 61 x 31 cm,
 - Length in extended position: 175 cm,
 - Overall: 155 x 76 x 138/168 cm,
- Maximum lifting capacity: 180 KG,
- Weight: 94 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

