



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ ARCO MEDI-MATIC 114.495 Gynecological chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **New upholstery with medical certificate in gray color,**
- **ARCO gynecological chair** is designed for gynecological examinations and procedures,
- **Technical specifications:**
 - Electric tilt adjustment of the chair within 50° by foot control,
 - Adjustable leg rests within: 240°,
 - Retractable waste bowl,
- **Included:**
 - Power cord,
 - Wired foot control,
 - Pull-out waste bowl with dimensions: 29 x 32 cm,
 - 2 pcs foot rests,
- Dimensions:
 - Headrest diameter: ø30 cm,
 - 2 segmented part of the backrest:
 - Upper: 56 x 30 cm,
 - Lower: 57 x 55 cm,
- Maximum load capacity: 180 KG,
- Weight: 115 KG,
- **Has a current inspection and is ready for use,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

