



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TOPCON CT-20D Non-contact tonometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 80 VA,
- **Menu in language: English,**
- **Technical specifications:**
 - Measurement range: 0 - 30 / 0 - 60 mmHg,
 - Adjustment range: -13 D to +14 D,
 - Working range: 11 mm,
 - Diagonal display: 5",
 - Built-in printer with paper width: 55 mm,
 - Manual canopy adjustment,
 - Manually adjustable chin,
- **Includes:**
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 30.5 x 48 x 45.5 cm,
- Weight: 20 KG,
- **Has a current inspection and is ready for operation,**
- **Issued Technical Passport (Service Report) valid for a period of 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

