



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS endoscope set (OEV261H, CV-180, CLV-180, MAJ-1430, WM-NP1 (Reconditioned))

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **The set includes:**
 - **OLYMPUS LCD medical monitor OEV261H:**
 - **Olympus OEV261H monitor** is a world-class monitor for displaying high-quality images through its performance,
 - The ability to display the image in the selected available shape preferred by the monitor operator,
 - The ability to customize your own preferred settings under you which allows for a more efficient procedure and reduced preparation time,
 - Made in Japan,
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - **Specifications:**
 - Screen: LCD 26",
 - Resolution: 1920 x 1200 px,
 - Contrast ratio: 1000:1,
 - Number of colors: 16.8 million,
 - Angle with acceptable viewing of the displayed image: 178°,
 - **Input/output connectors:** D-sub 15-pin×1) /VGA□640 × 480, fv = 60 Hz) / SVGA□800 × 600, fv = 60 Hz) / XGA□1024 × 768, fv = 60 Hz) / WXGA□1280 × 768, fv = 60 Hz) / SXGA□1280 × 1024, fv = 60 Hz) / UXGA□1600 × 1200, fv = 60 Hz) / WUXGA□1920 × 1200, fv = 60 Hz) / DVI-D connector,
 - **OLYMPUS CV-180 EVIS EXERA II image processor:**
 - Made in Japan,
 - Power supply 230 V,
 - Frequency 50-60 Hz,
 - Power 150 VA,
 - The unit has the ability to process high-definition HDTV images to provide the best image quality in endoscopy and laparoscopy,
 - Compatible with EVIS 100/130/140 series and EVIS EXERA 160 series endoscopes,
 - H output
 -
 - HD / SD SDI output for high quality video transfer,
 - Has a port for xD memory cards,
 - Foiled control panel,
 - **OLYMPUS CLV-180 EVIS EXERA II light source:**
 - 230 V power supply,
 - Frequency 50-60 Hz,
 - 500 VA power,
 - Adapted for flexible endoscopes,
 - Main lighting: xenon bulb MD 631 / 300 W,
 - Emergency lighting: halogen bulb 12 V / 35 W,
 - Foiled control panel,
 - Built-in air pump with electronic power control,
 - **PIGTALL OLYMPUS MAJ-1430 cable:**
 - Japanese manufacture,
 - Provides a connection between the video processor and the endoscope,
 - The MAJ-1430 is designed for use with the OLYMPUS CV-180 and CV-190 video processor and the 180 and 190 series endoscopes,
 - Spiral, flexible cable,
 - **OLYMPUS WM-NP1 Endoscope Tower:**
 - **The Olympus WM-NP1** is a stable and sturdy cart suitable for carrying endoscopes and delicate medical equipment,
 - **It features:**
 - Upper top with dimensions: 47 x 38 cm,
 - Lower top with dimensions: 47.5 x 56 cm,
 - 2 shelves with dimensions: 40 x 54 cm with manual height adjustment (every 1.5 - 2 cm) with Allen screws,
 - Handles for easy movement of the cart,
 - Holder for endoscopes with manual height adjustment,
 - Monitor holder with height adjustment within 20 cm (6 points every 3.5 cm),
 - Electrical installation, 12 power sockets for the medical equipment carried,
 - **Product features:**
 - Power switch on the front of the cart,
 - Detachable rear panel of the cart for easy cable management,
 - Power cords aesthetically concealed in panels on the sides of the cart
 -
 - Power cables aesthetically hidden in panels on the sides of the cart,
 - Corners of the shelves protected by protective caps, and the edges over the wheels with rubber,
 - Mobile on four casters with Ø 12 cm, two of them with brakes,
 - Half with Velcro fasteners designed for 2 bottles,
 - **Included:**
 - Power cords,
 - OLYMPUS MAJ-845 keyboard,
 - OLYMPUS MAJ-901 endoscope rinse bottle,
 - OLYMPUS VEK0N57 AC adapter,
 - Overall dimensions: 70 x 90 x 190 cm,
 - Weight: 160 KG,
 - **It has a current inspection and is ready for operation,**
 - **Issued Technical Passport (Service Report) valid for a period of 12 months,**
 - **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
 - Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

