



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



YESON E-Series 22l Black LED Autoclave (New)

New,

Modern Yeson E-series autoclave (LED) with a capacity of 22 liters,

European class B in accordance with the EN 13060 standard,

The device is designed for use in: hospitals, medical facilities (in clinics), including medical and dental offices, as well as laboratories. The autoclave can also be used for non-medical and veterinary applications.

Features and advantages of the E-series autoclave (LED):

Stainless steel chamber,

Built-in thermal printer,

Safety valve,

LED display,

Bacteriological filter,

Easy operation,

4 ready to use programs and 2 test programs,

Sterilization of packaged and unwrapped items,

Option to adjust sterilization time,

Test programs - Bowie & Dick/Helix Test, Vacuum test (leak test),

Option to maintain temperature until the next cycle,

Technical data:

Medical device: class IIb,

Chamber capacity: 22 liters,

Internal tanks for distilled and used water,

Loading method: basket with 3 trays,

Power supply: 230V \pm 10% , 50/60 Hz,

Rated power: 1950 W,

Maximum working pressure: 0.21 ~ 0.24 MPa / 2.1 bar ~ 2.4 bar,

Maximum operating temperature: 134°C ~ 137°C,

Dimensions with door closed: 680 mm x 505 mm x 425 mm,

Chamber dimensions: 250 x 350 mm,

Tray dimensions: top: 19.5 x 39.5 cm | middle: 19.5 x 39.5 cm | bottom: 19.5 x 39.5 cm,

Certificate of compliance with Directive 93/42/EEC (MDD),

Weight: 60 kg,

Has a current Technical Passport issued,

Warranty: 24 months,

