



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **YESON Autoclave YS-03L-E Series (New)**

New,

Modern Yeson YS-03L series autoclave with a capacity of 3 liters,

European class B in accordance with the EN 13060 standard,

The device is designed for use in: hospitals, medical institutions (in clinics), including medical and dental offices, as well as laboratories. The autoclave can also be used for non-medical and veterinary applications.

Features and advantages of the YS-03L series autoclave:

Class B fractional vacuum,

Built-in printer and SD card port,

Built-in clean and waste water tank,

Water quality sensor,

Touch control panel,

Multilingual menu including Polish,

Lock against accidental opening,

4 automatic operation programs,

2 test programs,

Additional drying program,

Pre-heating function,

Password protection function,

Loading method: rack with two trays,

Technical data:

Medical device: class IIb,

Chamber capacity: 3 liters,

Equipped with 4.3" touch screen display,

Power supply: 230V  $\pm$  10% , 50/60 Hz,

Rated power: 2900 W,

Maximum working pressure: 0.24 MPa,

Maximum operating temperature: 137°C,

Dimensions with door closed: 595 mm x 287 mm x 248 mm,

Chamber dimensions: 175 mm long, 55 mm high, 280 mm deep,

Tray dimensions: 18 x 28 cm,

Certificate of compliance with Directive 93/42/EEC (MDD),

Weight: 24 kg,

Has a current Technical Passport issued,

Warranty: 24 months,

