



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



YESON YS-2.5L ultrasonic cleaner (New)

New,

Model YESON YS-2,5L is designed for ultrasonic cleaning of tools, mechanisms, any objects that require thorough cleaning without disturbing the cleaned surface,

Features and advantages of the washer:

The ultrasonic cleaner is composed of a steel tub, to which transducers are attached that make the device vibrate at a specific frequency of 40kHz,

The vibrations propagate in the washing solution in the form of ultrasonic waves, producing thousands of tiny bubbles which, striking the cleaned object, remove any dirt,

Zones of vacuum and overpressure acting on the surfaces, crevices and recesses of the cleaned objects, lead to rapid removal of contaminants, and achieve an excellent cleaning effect,

The device is equipped with a heater that allows the cleaning fluid to be heated, and a timer that allows the device to be turned off after a programmed cleaning time,

The washer is designed for ultrasonic cleaning of tools, mechanisms, any objects that require thorough cleaning without disturbing the cleaned surface,

The housing is made of high-quality plastic,

Technical data:

Dimensions of the tank: 245*150*75mm,

Dimensions of the device: 320 x 270 x 270mm,

Total capacity: 2.5 liters,

Supply voltage 230 V, 50 Hz,

Ultrasound frequency: 40 kHz,

Ultrasound power: 60W,

Heater power: 100W,

Temperature range: 20~60°C,

Time setting range: 1-30 minutes,

Maximum length of objects: 250 mm,

Weight: 4 kg,

Has a current Technical Passport issued,

Warranty 12 months,

