



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



YESON 6L ultrasonic cleaner (New)

New,

Model YESON 6L is designed for ultrasonic cleaning of tools, mechanisms, any objects that require thorough cleaning without disturbing the cleaned surface,

The ultrasonic cleaner is composed of a steel tub, to which transducers are attached that vibrate the device at a specific frequency of 40kHz,

Intended use:

Specially designed for dental clinics, treatment rooms, beauty salons and veterinary offices, this ultrasonic cleaner combines advanced technology with functionality,

Unique Design:

Equipped with a 6-liter flat tank, it can effectively wash more instruments while saving cleaning fluids,

Intuitive Operation:

Features a capacitive touch panel with 4 buttons, ensuring ease of use and excellent water resistance. Safety is ensured by placing the electronic boards outside the tank,

Advanced Features:

Equipped with a tri-color LED display and touch control panel with time (5-30 min) and temperature (50-60 °C) settings, supported by two ceramic heaters,

Multifunctional Basket:

Includes a 17 x 28.5 cm stainless steel basket that also serves as a drainer. It allows the placement of larger utensils and is approximately 33 cm in diagonal,

Performance and Safety:

Equipped with two sets of ultrasonic transducers (60W, 38mm diameter), ensuring even activity throughout the tank. In addition, the durable drain valve and sturdy plastic housing add to the durability and safety of the device,

Technical specifications:

Frequency: 40 kHz,

Temperature settings: 50, 55, 60 °C,

Waterproof and shockproof,

Tank dimensions: 330*200*100mm,

Device dimensions: 560 x 310 x 230mm,

Basket: internal dimensions: 17 x 28.5 cm,

Total capacity: 6 liters,

Supply voltage: 230 V, 50 Hz,

Ultrasound frequency: 40 kHz,

Ultrasound power: 120W,

Heater power: 220W,

Temperature range: 50~60°C,

Time setting range: 5-30 minutes,

Weight: 9 kg,

Accessories:

Steel basket with drainer function,

Hopper lid,

Two stainless steel containers (0.5L) for small tools,

Hose for drain valve,

Has a current Technical Passport issued,

Warranty 12 months,

