



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



REICHERT AT 555 Automatic non-contact tonometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 87 VA,
- Measurement range: 0 - 60 mmHg,
- **The AT555** uses a patented, hands-free auto-adjustment system that allows any staff member to take an accurate IOP measurement in seconds,
- It automatically detects which eye is being measured, takes the measurement and displays the results on the 5" LCD screen,
- Measurement data can be easily printed at the touch of a button,
- Fully automated operation means fast, error-free readings,
- Simple to use, patients simply lean against the instrument's forehead rest while the operator presses one button to activate the fully automatic alignment and IOP measurement process,
- **Features:**
 - No chin rest, joystick or lifting controls. Patients simply lean against the forehead support, see the fixation target - and moments later experience the gentlest air intake available in any NCT. The procedure is quick, quiet and comfortable, with no patient discomfort,
 - Alignment is completely automated,
 - Easy-to-use icon-based operating system,
 - Clearly displays all measurement data,
 - Internal printer simplifies record keeping,
 - Electronic data transfer capability,
 - Maintenance-free, no bulb replacement required
 - The left/right sliding headrest, allows easy positioning of the patient's head for quick and convenient measurements.
 - The AT555 nosepiece contains an array of bright red LEDs. These provide the patient with an obvious visual cue for quick fixation and positioning, thus reducing measurement time,
 - User-friendly operator interface,
 - An "eye" icon on the display shows the operator how the instrument's automatic alignment system is setting up for measurement,
- **Selectable menus in ENG, FR, GER, ESP, PT, IT,**
- **Included:**
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 43 x 24 x 36 CM,
- Weight: 13.6 KG,
- **Has an up-to-date inspection and is ready for operation,**
- **Issued Technical Passport (Service Report) valid for a period of 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

