



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MedLED CHROME MC7 PRO Headlamp (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: through 2 new batteries,
- Frequency: 50/60 Hz,
- **The medLED CHROM headlamp** is an excellent portable light option. With an emphasis on comfort, light intensity and battery life, the CHROME lamp is designed to withstand the longest surgical procedures,
- **The CHROME MC7 PRO** provides maximum wearable comfort, first-class light, and comes with an extremely powerful battery,
- This lamp can illuminate any field in the operating room in terms of both intensity and spot size. It is also quiet, convenient and wireless,
- **Advantages:**
 - Cordlessness - no wires or batteries on the belt means no hassle,
 - Sterility - designed to meet the highest sterility requirements,
 - Quiet-Cool Technology - proprietary technology is used to cool the LED instantly, avoiding potential unnecessary heat buildup,
 - Adjustable light spot,
- **Technical specifications:**
 - Light intensity: 275,000 lux,
 - Estimated lifetime of up to 100,000 hours,
 - Color temperature: 5000 - 6500 K,
 - Battery life: 3.5,
 - Battery charging time: 6.5 hours - directly via USB,
 - Change batteries with speed, ease and without turning off the lamp,
 - High-speed USB connector for charging from a computer or power outlet,
 - Integrated battery status indicators,
 - Smooth digital dimming control,
 - Spot quality: sharp, round beam of uniform light
 - Full beam angle adjustment, left-to-right rotation, increased spot size range,
 - Mains plug versions: EU / US / UK / AU,
- **Included:**
 - USB power cord,
 - Charging adapters included: EU / US / UK / AU,
 - Dedicated case for storing the lamp,
 - 2 pcs. - keys for the case,
 - 2 new lithium-ion batteries,
 - Charger for 4 pcs of batteries (powered via USB),
- Total weight of the instrumentation: 4 KG,
- **It has a current inspection and is ready for use,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,**
 - **The warranty period does not include bulbs and other light sources installed in the device,**
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

