



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## MEMMERT WB7 Water bath (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1200 W,
- **Technical specifications:**
  - Bath entirely made of high-grade acid-resistant steel,
  - Usable capacity: 7 liters,
  - Temperature control in the range: from ambient temperature to 95°C,
  - Adjustable timer in the range: from 1 minute to 99 hours,
  - Ability to set delayed start of heating,
  - Built-in digital display,
  - Simple and intuitive operation,
- Dimensions:
  - Overall: 49 x 38 x 24 cm,
  - Chamber: 25 x 24 x 15.5 cm,
- Weight: 11 KG,
- **Has up-to-date inspection and is ready for operation,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
  - **6-months for domestic market (Poland),**
  - **3-months for the international market,**
  - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

