



## **AESCLAP FLUID 3 PG-070 Laparoscopic irrigation pump (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 36 VA,
- **Fluid 3** is a high-performance irrigation pump for laparoscopy,
- **Technical specifications:**
  - Max. Flushing and suction power: 400mmHg,
  - Flow rate: 0- 1.8 ml / min,
  - Very easy to use system,
  - Pressure controlled pump,
- **Included:**
  - Power cord,
- Dimensions: 14 x 22 x 32 cm,
- Weight: 6 KG,
- **Has a current inspection and is ready for use,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
  - **6-months for domestic market (Poland),**
  - **3-months for the international market,**
  - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

