



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NOUVAG DISPENSER DP 20 Tumescient anesthesia pump (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 20 VA,
- **DISPENSER DP 20** is a universal pump used for tumescent anesthesia in endovascular, liposuction and other procedures,
- Intuitive to use,
- Equipped with a sturdy handle and a metal stand for hanging a bottle of tumescent fluid,
- Infinite, precise flow rate control with a potentiometer from 0 - 100%,
- Vario foot switch allows you to turn on, off and adjust the flow speed without interrupting your work,
- **Technical specifications:**
 - Flow rate: 0 - 22 L / h,
 - Maximum pressure: 2.5 bar,
- **Included:**
 - User manual in EN, DE, FR (PDF),
 - 35 cm high tripod for hanging the tumescent fluid bottle,
 - Vario foot switch with 250 cm long cable,
 - Power cord,
- Dimensions:
 - Pump with tripod: 26.5 x 28 x 50 cm,
 - Foot switch: 22 x 7.5 x 6 cm,
- Weight: 6 KG,
- **It has a current inspection and is ready for operation,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

