



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



THERMO SCIENTIFIC SORVALL ST8 Laboratory centrifuge (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the PRC, (2017),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 310 W,
- **SORVALL ST 8R** with cooling, is a small table-top centrifuge that finds application in both diagnostic and scientific laboratories. The wide range of available rotors and adapters allows centrifugation of vessels and tubes from 0.2 ml to 145 ml at speeds up to 16,000 rpm, and the unique Auto-lock system used allows quick and safe rotor changes without the need for tools,
- **Key features:**
 - Acceleration up to 24,328 x g,
 - Cooling: Yes,
 - Auto-lock quick rotor change system,
 - Rapid acceleration and deceleration of rotors,
 - Microprocessor control,
 - Ability to save 4 user programs,
 - Quick selection buttons dedicated to basic programs,
 - Easy-to-read digital display,
 - Indication of RPM speed or RCF acceleration on the display,
 - Automatic closing and opening of the lid,
 - Automatic rotor recognition,
 - Rotor unbalance protection system,
 - Operating time clock,
 - Pulse spinning capability,
 - Quiet operation during centrifugation,
 - Induction motor with overheating protection,
 - Soft start and soft stop rotor function,
 - 14 available rotors for tubes from 0.2ml to 145ml, MTP plates and hematocrit capillaries,
- **Specifications:**
 - Speed: 16,000 RPM (**no rotor included - can be retrofitted at additional cost**),
 - Acceleration: 30,279 x g
 - Running time: 10 s - 99 min, continuous mode
 - Temperature control: -10 °C to +40 °C
 - Programming: 99 programs
- **English language menu,**
- **Can be retrofitted with selected accessories (cups / adapters / test tubes) at additional cost),**
- **Included:**
 - User manual in English (PDF),
 - Power cord,
 - **No rotor included - can be retrofitted at additional cost,**
- Dimensions: 32 x 46 x 67 cm,
- Weight: 31 KG,
- **Has a current inspection and is ready for operation,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

