



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



THERMOGEAR CHILL BUSTER 8001 Patient warmer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: through 1 new built-in 12 V / 7 Ah battery,
- Frequency: 50/60 Hz,
- Power: 5 W,
- **Specifications:**
 - Infinitely adjustable temperature range: 20 - 40.5°C,
 - Has visual and audible alarms - POWER light, charging light, low battery light,
 - Maximum length of operation on a charged battery: up to 3 hours,
 - Can be hung on the bed frame,
- **Included:**
 - Heating mat with dimensions: 76 x 152 cm,
 - 2 pcs. - covers,
 - Power supply,
 - Power cable,
 - 12V car charger to power the device,
 - Portable dedicated carrying case,
 - User manual in English,
- Dimensions: 23 x 42 x 20 cm,
- Weight: 7 KG,
- **Has a current inspection and is ready for use,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

