



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARESTREAM RVG 6100 Radiovisiograph (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- High quality intraoral imaging,
- With its real resolution of 20 lp/mm, it guarantees high image quality. **RVG 6100** is a top quality product created for all those who expect the highest image precision. The system uses a sensor shape that allows easier and more precise positioning. Its rounded edges are more comfortable for the patient,
- The Carestream **RVG 6100** digital radiovisiography system allows the operator to take and analyze radiographs easily and quickly,
- Its reduced size facilitates proper and comfortable positioning of the sensor in children. The use of a unique technology with an optical fiber layer guarantees precise diagnosis, as well as unprecedented tolerance of exposure times. This large exposure time tolerance allows the operator to take perfectly exposed images every time,
- **Features:**
 - Unique optical fiber technology for high quality imaging,
 - Great exposure tolerance guaranteeing perfect pictures every time,
 - Extremely user-friendly and logical Carestream software,
 - Superior imaging and diagnostic quality consistent over time,
- **Technical specifications:**
 - Sensor matrix: 1200 x 1600 px,
 - Active dimension of the sensor surface: 22.2 x 29.6 mm,
 - Image resolution: 20 lp/mm,
 - **Technology:**
 - CMOS,
 - Scintillator,
 - Optical fiber,
- **Included:**
 - English user manual,
 - Power cord,
 - TARGUS USB 2.0 HUB,
 - USB cable,
 - Dedicated carrying case,
 - Mounting accessories,
 - Set of adapters for outlets,
- **Has a current inspection and is ready for use,**
- **Issued Technical Passport (Service Report) valid for a period of 12 months,**
- **Warranty:**
 - **6-months for domestic market (Poland),**
 - **3-months for the international market,**
 - **Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



