



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **INVENTIS HARP BASIC Audiometer (Ex-Display)**

- New (ex-display),
- Real photos of the product,
- Made in Italy, (2019),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- HARP Basic diagnostic audiometer from Inventis - an audiometer that allows testing for air and bone conduction, with masking, with suprathreshold tests (Fowler, SISI, Stenger, Decay Ton, DLI, MHA),
- It has the ability to perform speech audiometry (by connecting a CD player with a disc of verbal material) and field testing (by connecting an active loudspeaker column),
- It is equipped with a large, easy-to-read color display with a menu in Polish, at the same time it is lightweight and very easy to use,
- Technical specifications:
  - Intensity:
    - AC: From -10 to 120 dB HL,
    - BC: From -10 to 90 dB HL,
    - Frequency range: 125 to 8000 Hz,
    - Masking: NBN / WN / SN,
  - Available tests:
    - Tonal audiometry,
    - Automatic threshold (modified Hughson-Westlake),
    - Speech audiometry (2 channels),
    - ALB,
    - SISI: automatic score - 1 dB increments (5 dB for familiarization),
    - Stenger, with pure tone or speech stimulation,
    - Tone decay with a duration of 60 or 120 seconds,
    - DLI, with 0 to 5 dB increments,
    - 2 independent Master Hearing Aid channels,
  - TFT LCD color display: 480 x 272 px,
  - Computer interface:
    - Connection: USB (no driver),
    - Compatible software: Inventis Maestro,
  - Database:
    - Up to 100 patients, AC and BC thresholds (tonal studies only),
- Included:
  - Power cord,
  - INVENTIS earphones,
  - SENNHEISER earphones,
  - Bone headphones,
  - Patient Button,
  - Microphone,
  - USB cable,
- Device dimensions: 32 x 32 x 9 cm,
- Weight of the device: 1.8 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

