



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FARUM A40 Oxygen cylinder regulator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Reducer is used to reduce the pressure of oxygen drawn from the cylinder,
- Integrated with a dispenser and humidifier designed for passive oxygen therapy - it is used to administer humidified oxygen or oxygen-air mixture to the patient,
- Supply pressure: 2 - 20 MPa,
- Flow control range: 0 - 25 l/min,
- Cylinder connection - 3/4" male thread,
- Dimensions: 12.5 x 10 x 9 cm,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

