



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET ALPHASTAR PRO 1132.21B Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 17 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- Mobile universal operating table with electrohydraulic drive and built-in emergency control panel on the table column,
- New upholstery (memory visco) with medical approval in black,
- Adjustable in multiple planes via wired remote control or via the panel on the column:
  - ON/OFF,
  - LOCK/UNLOCK - lock and unlock the table,
  - NORMAL/REV - orientation of the patient on the operating table,
  - TREND/REV - change of the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg positions in the range:  $\pm 30^\circ$ ,
  - The table is equipped with an auxiliary drive controlled by a hand-held remote control possible forward and backward movement,,
  - UP/DOWN - operating surface height (up/down) in range: 68.5 - 112 cm,
  - TILT L/TILT R - lateral tilt in the range of  $20^\circ$ ,
  - BACK UP/BACK DN - back segment position adjustment:  $-40^\circ / 0^\circ / +70^\circ$ ,
  - LEG UP/LEG DN - adjustment of leg segment: up/down,
  - FLEX/REFLEX - bending the top in half both ways: roof position / V position (Beach chair),
  - LEVEL - resetting: return to horizontal recliner position,
- Manual adjustment:
  - Headrest,
  - Two-segment leg rest for vertical tilt and side-to-side adjustment,
- Base housing made of durable plastic,
- The table is X - RAY translucent,
- Equipped with rails for attaching additional accessories,
- Table adapted for C-arm,
- Mobile, on double casters
- Manually removable headrest and foot segment,
- Includes:
  - Power cord,
  - Wired remote control,
- Possibility to retrofit the table at an additional cost with other accessories from our offer, examples shown in the last photos,
- Dimensions:
  - Headrest: 50 x 26 x 7 cm,
  - Backrest: 55 x 54 x 7 cm,
  - Seat: 55 x 42 x 7 cm,
  - Main leg rest segment: 55 x 52 x 7 cm,
  - Lower leg rest segment: 54 x 26 x 7 cm,
- Maximum load capacity: 454 KG,
- Weight: 300 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

