



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ELIGHT IPL laser therapy unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the PRC,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Application:
 - E-LIGHT:
 - Depilation,
 - Skin rejuvenation,
 - Pigmentation therapy,
 - Vascular therapy,
 - Acne therapy,
 - Wrinkle removal,
 - Laser:
 - Eyebrow tattoo removal,
 - Lip wrinkle removal,
 - Removal of wrinkles under the eye,
 - Tattoo removal,
 - Skin spot removal,
 - Removal of birthmarks on the skin,
 - IPL:
 - Depilation,
 - Skin rejuvenation,
 - Pigmentation therapy,
 - Vascular therapy,
 - Acne therapy,
 - Wrinkle removal,
 - RF:
 - RF Power: Up to 50 W,
 - PULSE / PWM,
 - Treatment time up to: 99 minutes,
 - 5-step head cooling control,
 - Modes of operation:
 - Facial areas,
 - Whole body (without face),
- Selectable menu in ENG, CHIN, FR, IT, ESP, RUS, ,
- Included:
 - Power cord,
 - 2 pcs heads,
- Dimensions: 51 x 70 x 151 cm,
- Weight: 45 KG,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

