



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



SANHE BEAUTY S&T SQ-2 Laser therapy unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the PRC, (2015),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1100 W,
- Selectable menu in language: ENG,
- Technical specifications:
 - Wavelength range:
 - 1064 nm,
 - 532 nm,
 - Energy level: 100 - 1500 mJ,
 - Pulse duration: 10 ± 3 ns,
 - Frequency: 1 - 10 Hz,
 - Speckle size can be adjusted by lens with pressure adjustment from 1 - 8 mm,
 - Spring arm with articulated joint with free adjustment,
 - Water cooling system of the machine for continuous use,
 - Operation and parameter adjustment via a 10.4" TFT color touch screen,
- Application:
 - Acne treatment,
 - Facelift,
 - Pigment removal,
 - Skin rejuvenation,
 - Tattoo removal,
 - Wrinkle removal,
- Included:
 - Power cord,
- Dimensions: 40 x 68 x 82 cm,
- Weight: 75 KG,
- Comes with a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

