



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EUROMI NIL SP6 Liposuction (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Belgium,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Aesthetic indications:
 - Figure shaping,
 - Liposuction of the buttock area,
 - Facial modeling, including restoring facial volume and shape, filling and smoothing some wrinkles, improving facial harmony,
 - Modeling, breast augmentation (e.g. increasing bust size, potential alternative to breast implants),
- Characterization:
 - Indication of the amount of infiltrated fluid per area,
 - Status indicators to guide the user through the infiltration or liposuction process,
 - Four modes of operation: infiltration/suction and pedal/continuous,
 - Infiltration cartridge connected to infiltration tubing set,
 - Compact control panel,
 - Easy to transport,
 - Improved working environment,
 - Noise reduction,
 - Overpressure control,
- Standard equipment:
 - Aspiration system (professional suction pump),
 - Infiltration system (quantitative infiltration pump),
 - SAFETY System,
 - 1x EVAMATIC 4 head (reconditioning) - can be exchanged for a new EVAMATIC 5 head at an extra charge,
 - 2x Foot switch,
 - 5x Pressure Hose,
 - 5x Suction hose,
 - 5x Connecting hose for infiltration,
 - 5x BAXTERS infiltration syringes (set),
 - 5x Medivac Flex 3L suction cartridge,
 - 2x Jar,
 - 1x Electric power cord,
 - 1x Cannula adjustment wrench,
 - 1x Cleaning wrench,
 - 1x Lipomatic® oil (100ml),
 - Cannula:
 - 1x Mercedes 4x300,
 - 1x 3G Accelerator 3x150,
 - 1x 3G Accelerator 3x250,
 - 1x Viterbo 4x250,
 - 1x Infiltration 2.5x200,
 - 1x Infiltration 2.5x300,
 - Additional equipment (at an extra charge):
 - Compressor DK50 2V Silent S DRY Monsun - in a silencing enclosure and with a dehumidifier,
 - Compressor DK50 2V DRY Monsun - without housing with dehumidifier,
 - MEDELA fat transfer kit (Fat Collector / Fat Transfer)
 - EUROMI NLF1 Fat Transfer Kit (Fat Collector / Fat Transfer)
 - EUROMI NLF2 Fat Transfer Kit (Fat Collector / Fat Transfer)
 - Consumables,
- Technical data of the suction machine:
 - JP 90:
 - Voltage (V~)/ Frequency: 200-240,
 - Maximum suction power (mbar): -850,
 - Power (W): 250,
 - Noise level db(A): 57,
 - KP 90:
 - Voltage (V~)/ Frequency: 100-115,
 - Maximum suction power (mbar): -700,
 - Power (W): 300,
 - Noise level db(A): 56,
- Front panel includes:
 - Control panel for quantity infiltration,
 - Air tubing sockets,
 - Infiltration syringes,
 - MEDI-VAC jars and cartridges,
 - Infiltration cartridges,
 - Pedicle sockets,
 - Silicone tubing socket used to connect the MEDI-VAC jar to the aspirator,
 - Emergency stop button,
 - Control panel:
 - infiltration / suction,
 - pedal / continuous mode,
 - Technical specifications:
 - Infrasonic drive,
 - Power supply:
 - 4A for 200 / 240V power supply,
 - 8A for 100 / 115V~ power supply,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3 months for international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

