



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TRIWINGS TWPDT Aesthetic phototherapy lamp (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in France,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 600 VA,
- TriWings® phototherapy represents the world's most powerful and effective LED photomodulation treatment system, combining multiple colors of light to induce interactions in tissues to reduce oxidative stress, reduce pain, reduce inflammatory reactions, promote healing of diseased or damaged skin, and inhibit and prevent the effects of aging processes. Phototherapy accelerates the healing and regeneration processes of tissues, significantly reducing the recovery time after procedures, including plastic surgery, aesthetic medicine, as well as in damage to the skin and subcutaneous tissue and in the course of a wide variety of diseases,
- LED light penetrates the skin to a depth of 8 mm, causing optimization of biochemical processes at the mitochondrial level. The radiation reaches the inner membrane, where cytochrome oxidase is stimulated into action. Activation of the cell's respiratory chain occurs, which increases cellular metabolism of melanocytes, keratinocytes, fibroblasts, and reduces free radicals (reduces oxidative stress) - such wide-ranging effects mean that LED phototherapy finds numerous applications in aesthetic medicine, plastic surgery and dermatology, among others,
- The LED system uses a spectrum of visible radiation from 380 to 780 nm and infrared 850 nm,
- Therapy indications:
 - Anti-aging skin treatments and PDT photostimulation,
 - Treatment of telangiectasias or pregnancy stretch marks,
 - Stimulation of the scalp and hair follicles to enhance hair growth and improve hair structure,
 - Baldness,
 - Scar treatment and pain therapy,
 - Topical therapy for acne, psoriasis, acquired vitiligo and chloasma,
 - Analgesic effects, muscular pains,
 - Eczema, allergies,
 - Burns,
 - Rash,
 - Juvenile acne, rosacea,
 - Scars, stretch marks,
 - Hard-to-heal wounds, bedsores,
 - Contusions, sports injuries,
 - Degenerative changes,
 - Supportive treatment of depression, stress and relaxing effect,
 - Complementing other therapies,
- Selectable menu in language: FR, ENG, PL, GER, ESP, IT,
- Included:
 - Power cord,
- Overall dimensions: 44 x 70 x 144 cm,
- Weight: 40 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
 - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

