



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET 1002.72A0-04 Headrest with MAQUET 1002.64A0 articulating handle (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photos of the product,
- Made in Germany,
- MAQUET 1002.72A0-04 is a two-segment, split headrest,
- It is used to support the patient's head in both supine and prone positions,
- Manual adjustment of the headrest in different positions,
- Adjustable width (inside) of the headrest in the range: 9 - 18 cm, each side separately,
- The 2 sections of the headrest are adjustable independently of each other,
- Suitable for use with MAQUET surgical table: 1130 / 1131 / 1132 / 1140 / 1150,
- Overall dimensions: 45 x 25 x 34 cm
- Weight: 7 KG,
- The headrest has normal signs of use that do not affect the comfort of use,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

