



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET 5100 Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- New upholstery (memory visco) with medical certificate in black,
- Technical specifications:
- Manual adjustment:
  - Height in the range: 77 - 102 cm,
  - Possibility to vertically raise the back segment in the range: 36 cm,
  - Tilt of the back segment: 90°,
  - Tilt of the upper back segment: 0 / -32°,
  - Tilt of upper seat segment: -70° / 0° / 32°,
  - Free side-to-side adjustment of the foot segment within the range: 90°,
  - Foot segment tilt -90° / 0° / 38°,
  - Trendelenburg position: 48°,
  - Anti-Trendelenburg position: 42°,
  - Side-to-side tilt: 32°,
- Dimensions:
  - Headrest: 21 x 24 cm,
  - 2-segment backrest: 48 x 52 cm,
  - 2 -segment seat: 48 x 35 cm,
  - 2 -segment foot rest:
    - Upper: 22 x 19 cm,
    - Lower: 22 x 56 cm,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

