



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **STORZ 250 TWIN light source 20113320 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 330 VA,  
Storz 250 TWIN  
is a halogen light source suitable for flexible endoscopes,  
Technical specifications:  
Halogen light,  
2 24 V 250 W incandescent lamps,  
Color temperature: 3000 K,  
Bulb life: 50 h,  
3-stage adjustment of light intensity with a knob,  
Bulbs protected by additional cover,  
Manual selection of the bulb to work with the included allen key hidden under the cover of the device,  
Input for the other end of the fiber optic cable protects against damage during a break in operation,  
Screwable fiber optic cable input allows its installation of an adapter corresponding to a fiber optic cable of another manufacturer  
(adapter not included)  
,  
Foiled control panel,  
Included:  
Power cord,  
Dimensions: 30.5 x 11 x 24.5 cm,  
Weight: 7.1 KG,  
Cat 02,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,  
The warranty period does not include bulbs and other light sources installed in the device,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

