



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OTV-SI Video image processor and light source (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Main and emergency lighting: 2 halogen bulbs 15 V / 150 W,
- Adapted for flexible and rigid endoscopes,
- Has the following functions:
 - Manual and automatic modes,
 - Contour enhancement for easy observation of tissues,
 - White balance,
 - Gradual adjustment of light intensity,
- It has the following inputs and outputs:
 - 2 x Y/C,
 - VIDEO OUT,
 - SYSTEM,
 - H. COPY,
 - VTR REMOTE,
 - Socket for keyboard connection,
- Keyboard on the control panel in language: English,
- Included:
 - OLYMPUS Keyboard Model MAJ-1124,
 - Power cord,
 - User manual in English (PDF),
- Dimensions: 43 x 30 x 17 cm,
- Weight: 13 KG,
- Cat 04,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
 - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

