



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **MEDIFA-HESSE 400308 Gynecological chair with DR. MACH 120F (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Technical specifications:
  - Free back tilt adjustment with the ability to unfold the chair to the reclining position,
  - Free tilt adjustment of the chair,
  - Free adjustable knee bolsters on both sides,
  - Built-in paper holder,
  - Freely extendable waste bowl,
  - Retractable patient step,
- Chair retrofitted with DR.MACH LED 120F wall mount lamp:
  - Technical specifications:
    - Lamp power: 40 W,
    - Light intensity: approximately 40,000 Lux,
    - Diameter of the illuminated field: 120 mm,
- Upholstery with medical approval in green color,
- Included:
  - Power cord,
  - 2x chair bolsters,
  - Wired foot control,
  - Waste bowl,
- Dimensions:
  - Back segment: 87 x 57 x 8 cm,
  - Seat: 55 x 32 x 8 cm,
  - Total length in reclining position: 140 cm,
- Weight:
  - Armchair: 135 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

