



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **BIOPHOTON DIGILED LED lamp (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Technical specifications:
  - 6 power levels:
    - 1st level - 0.5 W/cm<sup>2</sup>,
    - 2 level - 1 W/cm<sup>2</sup>,
    - 3 level - 1.5 W/cm<sup>2</sup>,
    - 4 level - 2 W/cm<sup>2</sup>,
    - 5th level - 2.5 W/cm<sup>2</sup>,
    - 6th level - 3 W/cm<sup>2</sup>,
  - 5 main modes:
    - Moisturizing,
    - Vaginal firming,
    - Firming the vulva,
    - Brightening of the intimate area,
    - Urinary incontinence,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

