



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MESOREGEN PLASMA GEL MAKER Hybrid filler manufacturing machine (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Realistic pictures of the product,
- (2019 r.),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 45 W,
- Thanks to the developed PPP filling technology, plasma taken from the blood is converted into a gel that is injected into the body,
- MesoREGEN is one of the indispensable and important equipment for the transformation of Blood Plasma into autologous collagen, which replaced very and adjusts the advanced heating and cooling system. Good features and no toxicity without side effects, it can also extract platelet-rich plasma (PRP) which is injected into the dermis, promotes skin renewal, promotes wrinkle elimination, improves facial appearance, diminishes pores and makes the skin soft again,
- Application:
 - Reduction of wrinkles and fine lines,
 - Improving skin elasticity and firmness,
 - Treatment of scars and stretch marks,
 - Improving skin color and reducing hyperpigmentation,
 - Treatment of baldness and improvement of scalp condition,
 - Skin recovery after sunbathing,
- Features:
 - Improving skin texture and elasticity,
 - Smoothing the skin and reducing the appearance of fine wrinkles,
 - Reduction of fine lines,
 - Improving skin firmness and tone,
 - Lightening of hyperpigmentation and pigmentation spots,
 - Reducing the visibility of surgical scars and acne scars,
 - Supporting skin regeneration processes,
 - Stimulating hair growth and improving the condition of the hair.
 - Stimulation of hair growth and improvement of scalp condition for trichology,
 - Skin reconstruction after excessive tanning,
 - Increasing skin hydration,
 - Improving skin density and volume,
 - Reducing the appearance of tired dark circles under the eyes,
 - Improving the overall appearance and healthy glow of the skin,
- Technical Specifications:
 - Display: Digital display,
 - Temperature range: RT~+90°C,
 - Capacity: 10 syringes,
- Included:
 - Power cord,
- Dimensions: 25 x 25 x 14 cm,
- Weight: 5.2 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

