



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **BENCHMARK MYTEMP H2200-H Mini incubator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the PRC,
- Power supply: 230 V / 12V,
- Frequency: 50/60 Hz,
- MYTEMP MINI H2200-H is a mini incubator with heating function,
- Basic features:
  - Incubator with heating,
  - Digital temperature control,
  - LCD display with current temperature,
  - For bottles up to 2000 ml,
- Specifications:
  - Temperature range: +5°C above ambient temperature to + 60°C,
  - Digital temperature control with increments of 1°C,
  - Temperature stability at 37 °C:  $\pm 0,5^{\circ}\text{C}$ ,
  - Temperature uniformity at 37°C:  $\pm 1^{\circ}\text{C}$ ,
  - Capacity: 20 liters,
- Included:
  - 230V power cord (supplying the internal outlet of the incubator),
  - Power supply + 12 V power cord,
  - Power cord (car 12 V),
  - 2 shelves with dimensions: 20.5 x 26 cm,
- Dimensions:
  - External: 33.5 x 37 x 48 cm,
  - Chamber: 26 x 23.5 x 33 cm,
- Weight: 7 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

