



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OE261H Medical monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Olympus OE261H monitor is a world-class monitor for displaying high quality images through its performance,
- The ability to display the image in the selected available shape preferred by the monitor operator,
- The ability to customize your own preferred settings under the operator which allows for a more efficient procedure and reduced preparation time,
- Technical specifications:
 - Screen: LCD 26",
 - Resolution: 1920 x 1200 px,
 - Contrast ratio: 1000:1,
 - Number of colors: 16.8 million,
 - Angle with acceptable viewing of the displayed image: 178°,
 - Input/output connectors:
 - D-sub 15-pin (x 1),
 - VGA (640 x 480, fv = 60 Hz),
 - SVGA (800 x 600, fv = 60 Hz),
 - XGA (1024 x 768, fv = 60 Hz),
 - WXGA (1280 x 768, fv = 60 Hz),
 - SXGA (1280 x 1024, fv = 60 Hz),
 - UXGA (1600 x 1200, fv = 60 Hz),
 - WUXGA (1920 x 1200, fv = 60 Hz),
 - DVI-D connector,
 - VGA (640 x 480, fv = 60 Hz),
 - SVGA (800 x 600, fv = 60 Hz),
 - XGA (1024 x 768, fv = 60 Hz),
 - WXGA (1280 x 768, fv = 60 Hz),
 - SXGA (1280 x 1024, fv = 60 Hz),
 - UXGA (1600 x 1200, fv = 60 Hz),
 - WUXGA (1920 x 1200, fv = 60 Hz),
- Selectable menus in ENG, GER, FR, CHIN, IT, ESP,
- The display has visible scratches on the anti-reflective coating that do not affect the comfort of use,
- Included:
 - User manual in English (PDF),
 - Power cord,
 - Monitor power adapter,
- Dimensions:
 - Monitor: 60 x 10 x 41 cm,
 - Power supply: 24 x 17 x 5 cm,
- Weight: 10.6 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

