



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **F1 FIMET OY Ophthalmic Chair (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Finland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- New upholstery with medical certificate in blue color (upholstery can be exchanged for another upholstery from our color chart at additional cost),
- Control by foot remote control:
  - Height adjustment in the range: from 48 to 94 cm,
  - Backrest adjustment in the range: from 0 to + 75°,
  - Leg rest adjustment: from -50 to 0°,
- Features:
  - Automatic positioning of the leg support along with the positioning movement of the backrest,
  - Program with setting the last position of the chair,
  - Possibility of manual removal of armrests from the chair,
  - Possibility of retracting the headrest,
  - Manual adjustment of the headrest within 180°,
  - Manual backrest extension within 22 cm,
  - Manual extension of the headrest within: 25cm,
  - Manual rotation of the chair right - left within 100° from the center position,
  - Possibility of manual tilt of the chair within the range of -3 to + 17°,
  - Possibility of manual tilt of the footrest in the range of -42 to - 30°,
  - Armrests with lateral recline function,
- Included:
  - User manual in English,
  - Power cord,
  - Power supply,
  - Foot switch,
- Dimensions:
  - Back segment: 42 x 64cm,
  - Seat: 62 x 55cm,
  - Leg rest: 50 x 52cm,
  - Total length of unfolded chair in straight position with headrest extended: 182cm,
  - Total width: 62cm,
  - Base:
    - Width: 51cm,
    - Length: 64cm,
- Maximum load capacity: 180kg,
- Weight: 118 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

