



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HUMPHREY 597 autorefractometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Technical specifications:
Measurement range:
VD of patient: 0 / 10.5 / 12 / 13.5 / 16.5 mm.,
Sphere (SPH): -17 D. to +20 D.,
Cylinder (CYL): From -7 D. to +7 D.,
Axis: Adjustable from 0° to 180°,
5" monitor,
Optotypes for children and adults,
Features built-in thermal printer,
Manual chin adjustment,
Electronic adjustment of the canopy and change of the examined side,
Output port: 2 x RS 232,
VIDEO output,
PCMCIA card slot and program card,
Selectable language menu: GER, ENG, FR, IT, ESP,
Included:
Power cord,
Dimensions: 32 x 51 x 48 cm,
Weight: 22.5 KG,
Has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

