



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BIOFRIDGE Portable Freezer BF-35-MIL-UPSD (Reconditioned)

Reconditioned (used),
Technical condition: good,
Visual condition: good,
Real photos of the product,
Manufactured by U.S.A., (2017),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 40 W,
BIOFRIDGE

is a 35-liter portable refrigerator/freezer is a compact model commonly used as a refrigerator for vaccines and blood transport, it is equipped with an internal fan, sturdy handles, a locking latch and a programmable digital thermostat,

Technical specifications:

Capacity: 35 l,
Temperature range in refrigerator mode: 2 °C to 8 °C,
Temperature range in freezer mode: 0 °C to -16 °C,
Built-in handles on both sides for free transportation,
Freezer flap manually lockable with latch,
Robust and compact design,

Included:

Power cord,
Car charger,
Internal basket,

The device has visible signs of use that do not interfere with free use,

Dimensions:

Internal: 35.5 x 28 x 38 cm,
External: 80 x 40 x 50 cm,
Inner basket: 31 x 23 x 34 cm,
Weight: 27 KG,

It has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),
3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

