



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Videoprinter SONY UP-895MD (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
Japanese production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Technical specification:  
Thermal printing,  
256 shades of gray,  
High print resolution: 1210 x 582 pixels,  
Multifunctional printing,  
Printing speed: 3.9 sec.,  
Intended for printing from CCTV industrial television systems and medical equipment, e.g. ultrasound,  
VIDEO input and output,  
Included:  
User manual in English (PDF),  
Power cable,  
Dimensions: 15.5x 27.5 x 10.5 cm,  
Weight: 3.3 KG,  
Cat 02,  
It has a current inspection and is ready to work,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,  
Financing options (Poland only): Installments, Leasing, Loan,

