



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Videoprinter SONY UP-895MD (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Technical specification:
Thermal printing,
256 shades of gray,
High print resolution: 1210 x 582 pixels,
Multifunctional printing,
Printing speed: 3.9 sec.,
Intended for printing from CCTV industrial television systems and medical equipment, e.g. ultrasound,
VIDEO input and output,
Included:
User manual in English (PDF),
Power cable,
Dimensions: 15.5x 27.5 x 10.5 cm,
Weight: 3.3 KG,
Cat 02,
It has a current inspection and is ready to work,
Technical Passport issued (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

