



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



OLYMPUS CLV-180 EVIS EXERA II light source (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Japanese production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 500 VA,
- **Technical specifications:**
 - Adapted to flexible endoscopes,
 - Main lighting: MD 631 / 300 W xenon bulb,
 - Emergency lighting: 12V/35W halogen bulb,
 - Foiled control panel,
 - Built-in air pump with electronic power control,
- **Included in the set:**
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 54 x 39 x 16 cm,
- Weight: 15 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

