



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



OLYMPUS CLV-180 light source (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 500 VA,
Technical specifications:
Adapted for flexible endoscopes,
Main lighting: xenon bulb MD 631 / 300 W,
Emergency lighting: 12 V / 35 W halogen bulb,
Foiled control panel,
Built-in air pump with electronic power control,
Included:
User manual in English (PDF),
Power cord,
Dimensions: 54 x 39 x 16 cm,
Weight: 15 KG,
Has an up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

